

Anti-Polr2D Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O15514

Item Number: CM0007287



Introduction

High-quality rabbit monoclonal antibody against human POLR2D (RPB4), validated for WB, IHC-P, ELISA. Excellent tool for studying RNA polymerase II transcription machinery. Research-use-only product from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	POLR2D Rabbit mAb
Remarks/Alias	RBP4; RPB4; RPB16; HSRBP4; HSRPB4; POLR2D
Species	Human
GeneID (Human)	5433
GeneID	5433
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-142 of human POLR2D (NP_004796.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15514
Protein Weight	16 kDa
Shipping	Ice bag