


Anti-Polr2F Polyclonal Antibody For Wb, Elisa - P61218

Item Number: CM0026636



Introduction

Rabbit polyclonal antibody targeting POLR2F (RPABC2), a common subunit of RNA polymerases I, II, III. Validated for WB and ELISA, reactive with human and rat samples. Ideal for transcription mechanism studies.

[Learn More](#)

Parameter	Value
Product Name	POLR2F Rabbit pAb
Remarks/Alias	RPB6; POLRF; RPC15; RPABC2; RPB14.4; HRBP14.4; RPABC14.4; POLR2F
Species	Human
GeneID (Human)	5435
GeneID	5435
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-127 of human POLR2F (NP_068809.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P61218
Protein Weight	14kDa
Shipping	Ice bag