

Anti-Polr2G Rabbit Polyclonal Antibody For Wb, Elisa - P62487

Item Number: CM0027110



Introduction

Rabbit polyclonal antibody against human RNA polymerase II subunit RPB7 (POLR2G). Validated for Western blot (WB) and ELISA; cross-reacts with human, mouse, and rat. Suitable for transcription and RNA polymerase II studies.

[Learn More](#)

Parameter	Value
Product Name	POLR2G Rabbit pAb
Remarks/Alias	RPB7; RPB19; hRPB19; hsRPB7; POLR2G
Species	Human
Gene ID (Human)	5436
Gene ID	5436
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human POLR2G (NP_002687.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62487
Protein Weight	19kDa
Shipping	Ice bag