


Anti-Polr2H Polyclonal Antibody For Wb, Ip, Elisa - P52434

Item Number: CM0026295



Introduction

Rabbit polyclonal antibody targeting POLR2H (RPB8/RPABC3), validated for WB, IP, and ELISA. Reacts with human, mouse, and rat samples. Ideal for RNA polymerase research and transcriptional studies.

[Learn More](#)

Parameter	Details
Product Name	POLR2H Rabbit pAb
Remarks/Alias	RPB8; RPB17; RPABC3; POLR2H
Species	Human
Gene ID	5437
Immunogen	Recombinant protein / Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human POLR2H (NP_006223.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52434
Protein Weight	17kDa
Shipping	Ice bag