


Anti-Polr2H Polyclonal Antibody For Wb, Ip, Elisa - P52434

Item Number: CM0015474



Introduction

Rabbit polyclonal antibody against DNA-directed RNA polymerases I, II, and III subunit RPABC3 (POLR2H). Suitable for Western Blot, Immunoprecipitation, and ELISA in Human and Mouse samples.

[Learn More](#)

Parameter	Value
Product Name	POLR2H Rabbit pAb
Remarks/Alias	RPB8; RPB17; RPABC3; POLR2H
Species	Human
GeneID (Human)	5437
GeneID	5437
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human POLR2H (NP_006223.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52434
Protein Weight	17kDa
Shipping	Ice bag