


Anti-Polr2J Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P52435

Item Number: CM0023481



Introduction

Rabbit polyclonal antibody against human POLR2J (RPB11-a), validated for WB, IF/ICC, ELISA. Suitable for RNA polymerase II transcription studies. Cross-reacts with mouse; nuclear localization and ubiquitous expression with high levels in heart and muscle.

[Learn More](#)

Parameter	Value
Product Name	POLR2J Rabbit pAb
Remarks/Alias	RPB11; RPB11A; RPB11m; hRPB14; POLR2J1; POLR2J
Species	Human
GeneID (Human)	5439
GeneID	5439
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-117 of human POLR2J (NP_006225.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	P52435
Protein Weight	13kDa
Shipping	Ice bag