


Anti-Polr2L Polyclonal Antibody For Wb, Ihc-P, Elisa - P62875

Item Number: CM0023763



Introduction

POLR2L Rabbit pAb, validated for WB, IHC-P, ELISA. Reacts with Human, Mouse, Rat. Targets RNA polymerase subunit RPABC5; ideal for transcription and chromatin studies.

[Learn More](#)

Parameter	Value
Product Name	POLR2L Rabbit pAb
Remarks/Alias	RBP10; RPB10; RPABC5; RPB7.6; hRPB7.6; RPB10beta; POLR2L
Species	Human
Gene ID (Human)	5441
Gene ID	5441
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-67 of human POLR2L (NP_066951.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62875
Protein Weight	8kDa
Shipping	Ice bag