

Anti-Polr3C Polyclonal Antibody For Wb, Elisa - Q9Bui4

Item Number: CM0017430



Introduction

POLR3C Rabbit pAb validated for WB and ELISA. Reacts with human, mouse, rat. Targets RNA Polymerase III subunit RPC3 (RPC62), a DNA sensor and transcription component. Ideal for innate immunity and transcription research.

[Learn More](#)

Parameter	Value
Product Name	POLR3C Rabbit pAb
Remarks/Alias	C82; RPC3; RPC62; POLR3C
Species	Human
Gene ID (Human)	10623
Gene ID	10623
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human POLR3C (NP_006459.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BUI4
Protein Weight	61kDa
Shipping	Ice bag