


Anti-Polr3F Polyclonal Antibody For Wb And Elisa - Q9H1D9

Item Number: CM0017944



Introduction

Rabbit polyclonal antibody targeting human POLR3F (RPC6, RPC39), validated for WB and ELISA. Cross-reacts with mouse and rat. Detects a 36kDa band, suitable for RNA polymerase III research and innate immunity studies.

[Learn More](#)

Parameter	Value
Product Name	POLR3F Rabbit pAb
Remarks/Alias	C34; RPC6; RPC39; IMD101; POLR3F
Species	Human
GeneID (Human)	10621
GeneID	10621
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human POLR3F (NP_006457.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H1D9
Protein Weight	36kDa
Shipping	Ice bag