


Anti-Polr3K Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y2Y1

Item Number: CM0028030



Introduction

Rabbit polyclonal antibody targeting human POLR3K, suitable for Western blot and ELISA. Detects a key RNA polymerase III subunit involved in transcription fidelity and innate immune sensing.

[Learn More](#)

Parameter	Value
Product Name	POLR3K Rabbit pAb
Remarks/Alias	C11, HLD21, My010, RPC10, RPC11, RPC12.5, C11-RNP3, POLR3K
Species	Human
Gene ID	51728
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-108 of human POLR3K (NP_057394.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9Y2Y1
Protein Weight	12kDa
Shipping	Ice bag