

Anti-Polrmt Rabbit Polyclonal Antibody For Wb, If/Icc, Ip, Elisa - O00411

Item Number: CM0022769



Introduction

Rabbit anti-POLRMT polyclonal antibody validated for WB, IF/ICC, IP, ELISA. Reacts with human and rat mitochondrial RNA polymerase, enabling mitochondrial transcription and replication studies.

[Learn More](#)

Parameter	Value
Product Name	POLRMT Rabbit pAb
Remarks/Alias	APOLMT; MTRNAP; MTRPOL; COXPD55; h-mtRPOL; POLRMT
Species	Human
Gene ID (Human)	5442
Gene ID	5442
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 42-230 of human POLRMT (NP_005026.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Rat
SWISS	O00411
Protein Weight	139kDa
Shipping	Ice bag