

Anti-Pomc Polyclonal Antibody For Wb, If/Icc, Elisa - P01189

Item Number: CM0027594



Introduction

Anti-POMC polyclonal antibody validated for Western blot, IF/ICC, and ELISA in human, mouse, and rat samples. Targets pro-opiomelanocortin, a precursor of pituitary hormones involved in energy balance and stress response.

[Learn More](#)

Parameter	Value
Product Name	POMC Rabbit pAb
Remarks/Alias	LPH; MSH; NPP; POC; ACTH; CLIP; OBAIRH; POMC
Species	Human
Gene ID	5443
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human POMC (NP_000930.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P01189
Protein Weight	29kDa
Shipping	Ice bag