


Anti-Pomgnt1 Rabbit Polyclonal Antibody For Wb, Elisa - Q8Wza1

Item Number: CM0030238



Introduction

High-quality rabbit polyclonal antibody against human POMGNT1 (GnT I.2), validated for western blot and ELISA. Reacts with human, mouse, and rat samples. Ideal for O-mannosyl glycosylation research and dystroglycanopathy studies.

[Learn More](#)

Parameter	Value
Product Name	POMGNT1 Rabbit pAb
Remarks/Alias	MEB; RP76; GNTI.2; LGMD20; GnT I.2; LGMDR15; MGAT1.2; gnT-I.2; POMGNT1
Species	Human
Gene ID (Human)	55624
Gene ID	55624
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 411-660 of human POMGNT1 (NP_060209.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WZA1
Protein Weight	75kDa
Shipping	Ice bag