

Anti-Pomgnt2 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nat1

Item Number: CM0017937



Introduction

Rabbit polyclonal antibody targeting human POMGNT2, validated for Western blot and ELISA. Cross-reacts with mouse and rat. Suitable for studying O-mannosylation, dystroglycan glycosylation, and congenital muscular dystrophy mechanisms.

[Learn More](#)

Parameter	Value
Product Name	POMGNT2 Rabbit pAb
Remarks/Alias	AGO61; GTDC2; MDDGA8; MDDGC8; C3orf39; POMGNT2
Species	Human
GeneID (human)	84892
GeneID	84892
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-340 of human POMGNT2 (NP_116195.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8NAT1
Protein Weight	67kDa
Shipping	Ice bag