

Anti-Pomt1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Y6A1

Item Number: CM0017771

Introduction

Rabbit polyclonal antibody against human POMT1 (Protein O-mannosyl-transferase 1). Validated for Western blot (WB) and ELISA, cross-reacts with human and mouse. Predicted band size ~85kDa. Ideal for dystroglycanopathy and glycosylation research.

[Learn More](#)

Parameter	Value
Product Name	POMT1 Rabbit pAb
Remarks/Alias	RT; LGMD2K; MDDGA1; MDDGB1; MDDGC1; LGMDR11; POMT1
Species	Human
Gene ID	10585
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 310-550 of human POMT1 (NP_001070833.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9Y6A1
Protein Weight	85kDa
Shipping	Ice bag