


Anti-Pon1 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P27169

Item Number: CM0009655



Introduction

Rabbit monoclonal antibody against human PON1 (Serum paraoxonase/arylesterase 1). Validated for WB, IF-P, IHC-P, ELISA. Cross-reacts with mouse and rat. Detects a 40 kDa protein (SWISS: P27169).

[Learn More](#)

Parameter	Value
Product Name	PON1 Rabbit mAb
Remarks/Alias	ESA; PON; MVCD5; PON1
Species	Human
GeneID (Human)	5444
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human PON1 (P27169).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27169
Protein Weight	40kDa
Shipping	Ice bag