


Anti-Pon1 Rabbit Polyclonal Antibody For Wb, Elisa - P27169

Item Number: CM0018182



Introduction

Rabbit polyclonal antibody against PON1 (P27169) for WB and ELISA. Cross-reacts with human, mouse, rat. Ideal for studying serum paraoxonase/arylesterase 1 function in organophosphate hydrolysis and HDL-associated antioxidant activity.

[Learn More](#)

Parameter	Value
Product Name	PON1 Rabbit pAb
Remarks/Alias	ESA; PON; MVCD5; PON1
Species	Human
Gene ID	5444
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 9-109 of human PON1 (NP_000437.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P27169
Protein Weight	40kDa
Shipping	Ice bag