

Anti-Pon2 Rabbit Polyclonal Antibody For Wb And Elisa - Q15165

Item Number: CM0022043

Introduction

Rabbit polyclonal antibody against human PON2 (serum paraoxonase/arylesterase 2), validated for WB and ELISA. Targets a 39 kDa membrane protein involved in antioxidant activity and lipid peroxidation prevention.

[Learn More](#)

Parameter	Value
Product Name	PON2 Rabbit pAb
Remarks/Alias	PON2
Species	Human
GenelD	5445
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human PON2 (NP_000296.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15165
Protein Weight	39kDa
Shipping	Ice bag