

Anti-Pon3 Polyclonal Antibody For Wb, Elisa - Q15166

Item Number: CM0019795

Introduction

High-quality rabbit polyclonal antibody targeting human PON3, validated for WB and ELISA applications. Ideal for investigating serum paraoxonase/lactonase 3 function and statin prodrug hydrolysis. Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	PON3 Rabbit pAb
Remarks/Alias	PON3
Species	Human
Gene ID (Human)	5446
Gene ID	5446
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 250-350 of human PON3 (NP_000931.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15166
Protein Weight	40kDa
Shipping	Ice bag