

Anti-Pop4 Polyclonal Antibody For Wb, Ihc-P, Elisa - 095707

Item Number: CM0027410



Introduction

Rabbit polyclonal antibody targeting POP4 (RPP29), a ribonuclease P subunit involved in tRNA maturation. Validated for WB, IHC-P, ELISA in human, mouse, rat. Ideal for tRNA processing research.

[Learn More](#)

Parameter	Value
Product Name	POP4 Rabbit pAb
Remarks/Alias	RPP29; POP4
Species	Human
GeneID (human)	10775
GeneID	10775
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human POP4 (NP_006618.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95707
Protein Weight	25kDa
Shipping	Ice bag