

Anti-Pop4 Rabbit Polyclonal Antibody For Wb, Elisa - 095707

Item Number: CM0016203



Introduction

High-quality anti-POP4 rabbit polyclonal antibody validated for WB and ELISA, targeting human and rat ribonuclease P subunit p29. Ideal for tRNA processing and nucleolar research.

[Learn More](#)

Parameter	Value
Product Name	POP4 Rabbit pAb
Remarks/Alias	RPP29; POP4
Species	Human
GeneID (Human)	10775
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-220 of human POP4 (NP_006618.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	O95707
Protein Weight	25kDa
Shipping	Ice bag