

Anti-Pop5 Rabbit Polyclonal Antibody For Wb And Elisa - Q969H6

Item Number: CM0021789



Introduction

High-quality rabbit polyclonal antibody targeting human POP5, a key subunit of ribonuclease P and MRP complexes involved in tRNA and rRNA maturation. Validated for WB and ELISA; cross-reactivity with human and rat. - Q969H6

[Learn More](#)

Parameter	Specification
Product Name	POP5 Rabbit pAb
Remarks/Alias	RPP2; RPP20; hPop5; HSPC004; POP5
Species	Human
Gene ID	51367
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human POP5 (NP_057002.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q969H6
Protein Weight	19kDa
Shipping	Ice bag