

Anti-Pop7 Rabbit Polyclonal Antibody For Wb, Elisa - 075817

Item Number: CM0024924



Introduction

Validate endogenous POP7 expression with this rabbit polyclonal antibody (075817) for Western blot and ELISA; suitable for investigating ribonuclease P/MRP complexes in tRNA and rRNA processing.

[Learn More](#)

Parameter	Value
Product Name	POP7 Rabbit pAb
Remarks/Alias	RPP2; RPP20; 0610037N12Rik; POP7
Species	Human
GeneID (Human)	10248
GeneID	10248
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human POP7 (NP_005828.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	075817
Protein Weight	16kDa
Shipping	Ice bag