

Anti-PorcN Rabbit Polyclonal Antibody For Wb And Elisa - Q9H237

Item Number: CM0015662



Introduction

High-quality rabbit polyclonal antibody targeting human PORCN, a Wnt signaling palmitoleoyltransferase, validated for WB and ELISA with cross-reactivity to mouse and rat. Ideal for Wnt pathway and disease research.

[Learn More](#)

Parameter	Value
Product Name	PORCN Rabbit pAb
Remarks/Alias	PPN; DHOF; FODH; MG61; PORC; PORCN
Species	Human
Gene ID (Human)	64840
Gene ID	64840
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 220-320 of human PORCN (NP_982299.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H237
Protein Weight	52kDa
Shipping	Ice bag