

Anti-Pot1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nux5

Item Number: CM0015814



Introduction

Rabbit polyclonal antibody against POT1, validated for WB and ELISA. Reacts with human and mouse POT1, a shelterin complex component essential for telomere length regulation and protection.

[Learn More](#)

Parameter	Value
Product Name	POT1 Rabbit pAb
Remarks/Alias	GLM9; CMM10; HPOT1; POT1
Species	Human
GeneID (Human)	25913
GeneID	25913
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 497-634 of human POT1 (NP_056265.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NUX5
Protein Weight	71kDa
Shipping	Ice bag