

Anti-Pou1F1 Rabbit Polyclonal Antibody For Wb And Elisa - P28069

Item Number: CM0024466



Introduction

High-quality rabbit polyclonal antibody against human POU1F1 (PIT1/GHF-1), validated for Western blot and ELISA. Ideal for pituitary development and transcription factor research.

[Learn More](#)

Parameter	Value
Product Name	POU1F1 Rabbit pAb
Remarks/Alias	PIT1; CPHD1; GHF-1; Pit-1; POU1F1a; POU1F1
Species	Human
GeneID (Human)	5449
GeneID	5449
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 124-273 of human POU1F1 (NP_000297.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P28069
Protein Weight	33kDa
Shipping	Ice bag