

Anti-Pou1F1 Rabbit Polyclonal Antibody For Wb, If-P, Ihc-P, Elisa - P28069

Item Number: CM0023812



Introduction

Rabbit polyclonal antibody against POU1F1 (PIT1/GHF-1), suitable for Western blot, immunofluorescence, immunohistochemistry, and ELISA. Detects human, mouse, and rat POU1F1, a pituitary-specific transcription factor involved in growth hormone and prolactin regulation. Ideal for pituitary development and endocrine research.

[Learn More](#)

Parameter	Specification
Product Name	POU1F1 Rabbit pAb
Remarks/Alias	PIT1; CPHD1; GHF-1; Pit-1; POU1F1a
Species	Human
Gene ID (Human)	5449
Gene ID	5449
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P28069
Protein Weight	33kDa
Shipping	Ice bag