

Anti-Pou2Af1 Polyclonal Antibody For Wb And Elisa - Q16633

Item Number: CM0029176



Introduction

Rabbit polyclonal antibody against POU2AF1, a B-cell-specific transcriptional coactivator essential for germinal center formation. Validated for WB and ELISA, cross-reacts with human, mouse, rat.

[Learn More](#)

Parameter	Value
Product Name	POU2AF1 Rabbit pAb
Remarks/Alias	BOB1; OBF1; OCAB; OBF-1; POU2AF1
Species	Human
Gene ID (Human)	5450
Gene ID	5450
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-256 of human POU2AF1 (NP_006226.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16633
Protein Weight	27kDa
Shipping	Ice bag