

Anti-Pou2F2 (Oct2) Rabbit Polyclonal Antibody For Wb And Elisa - P09086

Item Number: CM0021055



Introduction

Rabbit polyclonal antibody against human POU2F2 (OCT2, OTF2), validated for Western blot and ELISA. Cross-reacts with mouse and rat. Ideal for transcription factor studies in B-cell and immunoglobulin regulation research.

[Learn More](#)

Parameter	Value
Product Name	POU2F2 Rabbit pAb
Remarks/Alias	OCT2; OTF2; Oct-2; POU2F2
Species	Human
GeneID (human)	5452
GeneID	5452
Immunogen	Recombinant fusion protein containing a sequence within amino acids 400-479 of human POU2F2 (NP_001193954.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09086
Protein Weight	51kDa
Shipping	Ice bag