


Anti-Pou2F3 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uki9

Item Number: CM0013788



Introduction

Rabbit polyclonal antibody against human POU2F3 (Oct-11), suitable for WB and ELISA, cross-reacts with mouse. Targets a transcription factor involved in keratinocyte differentiation and tuft cell function. Ideal for dermatological and mucosal immunity research.

[Learn More](#)

Parameter	Value
Product Name	POU2F3 Rabbit pAb
Remarks/Alias	PLA1; OCT11; PLA-1; Epoc-1; OCT-11; OTF-11; Skn-1a; POU2F3
Species	Human
Gene ID (Human)	25833
Gene ID	25833
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human POU2F3 (NP_055167.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9UKI9
Protein Weight	47kDa
Shipping	Ice bag