

Anti-Pou3F2 Polyclonal Antibody For Wb, Elisa - P20265

Item Number: CM0023329



Introduction

Affinity-purified rabbit polyclonal antibody against human POU3F2 transcription factor (BRN2, OCT7). Validated for Western blot and ELISA; targets neuroectodermal lineage. Ideal for neuronal differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	POU3F2 Rabbit pAb
Remarks/Alias	BRN2; OCT7; OTF7; OTF-7; POUF3; brn-2; oct-7; N-Oct3; POU3F2
Species	Human
Gene ID (Human)	5454
Gene ID	5454
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-68 of human POU3F2 (NP_005595.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P20265
Protein Weight	47kDa
Shipping	Ice bag