

Anti-Pou3F4 Polyclonal Antibody For Wb And Elisa - P49335

Item Number: CM0019362

Introduction

Rabbit polyclonal antibody against human POU3F4 (BRN-4), validated for WB and ELISA. Detects human, mouse, and rat POU3F4. Suitable for neuroscience research and transcription factor studies.

[Learn More](#)

Parameter	Value
Product Name	POU3F4 Rabbit pAb
Remarks/Alias	BRN4; DFN3; OTF9; BRN-4; DFNX2; OCT-9; OTF-9; BRAIN-4; POU3F4
Species	Human
Gene ID (Human)	5456
Gene ID	5456
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human POU3F4 (NP_000298.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49335
Protein Weight	39kDa
Shipping	Ice bag