

Anti-Pou4F1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q01851

Item Number: CM0023253

Introduction

Rabbit polyclonal antibody against human POU4F1 (BRN3A), cross-reactive with mouse and rat. Validated in WB, IF/ICC, and ELISA. Ideal for neuronal differentiation, apoptosis, and transcription regulation studies.

[Learn More](#)

Parameter	Value
Product Name	POU4F1 Rabbit pAb
Remarks/Alias	BRN3A, RDC-1, ATITHS, Oct-T1, brn-3A
Species	Human
Gene ID (Human)	5457
Gene ID	5457
Immunogen	Recombinant protein; A synthetic peptide corresponding to amino acids 1-100 of human POU4F1 (NP_006228.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q01851
Protein Weight	43kDa
Shipping	Ice bag