

Anti-Pou4F3 Polyclonal Antibody For Wb, Elisa - Q15319

Item Number: CM0025378

Introduction

Rabbit polyclonal antibody targeting human POU4F3 (BRN3C), validated for Western blot and ELISA. Suitable for research on auditory system development and inner ear hair cell differentiation.

[Learn More](#)

Parameter	Value
Product Name	POU4F3 Rabbit pAb
Remarks/Alias	BRN3C; DFNA15; DFNA42; DFNA52; POU4F3
Species	Human
Gene ID	5459
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human POU4F3 (NP_002691.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q15319
Protein Weight	37kDa
Shipping	Ice bag