

Anti-Pou6F1 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q14863

Item Number: CM0028473



Introduction

Rabbit polyclonal antibody targeting human POU6F1 (BRN5, MPOU, TCFB1). Suitable for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience and transcription regulation research.

[Learn More](#)

Parameter	Value
Product Name	POU6F1 Rabbit pAb
Remarks/Alias	BRN5; MPOU; TCFB1; POU6F1
Species	Human
GeneID (Human)	5463
GeneID	5463
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-301 of human POU6F1 (NP_002693.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14863
Protein Weight	33kDa
Shipping	Ice bag