

Anti-Pp1 Beta Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62140

Item Number: CM0011260

Introduction

Recombinant rabbit monoclonal antibody targeting human PP1 beta catalytic subunit (P62140). Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Essential for cell division, glycogen metabolism, and MAPK pathway research.

[Learn More](#)

Parameter	Value
Product Name	PP1 beta Rabbit mAb
Remarks/Alias	MP; PP1B; PP1c; NSLH2; PP-1B; PPP1CD; PP1beta; PPP1beta; HEL-S-80p; PP1 beta
Species	Human
GeneID (human)	5500
GeneID	5500
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 228-327 of human PP1 beta (P62140).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62140
Protein Weight	37kDa
Shipping	Ice bag