

Anti-Pp1 Beta Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62140

Item Number: CM0013083

Introduction

Rabbit polyclonal antibody against human PP1 beta (PPP1CB) for WB, IHC-P, IF/ICC, ELISA applications. Cross-reacts with mouse and rat. Detects 37 kDa protein. Validated for cell cycle, glycogen metabolism, and circadian rhythm studies.

[Learn More](#)

Parameter	Value
Product Name	PP1 beta Rabbit pAb
Remarks/Alias	MP; PP1B; PP1c; NSLH2; PP-1B; PPP1CD; PP1beta; PPP1beta; HEL-S-80p; PP1 beta
Species	Human
GeneID	5500
Immunogen	Synthetic peptide (a synthetic peptide corresponding to a sequence within amino acids 227-327 of human PP1 beta (NP_002700.1))
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62140
Protein Weight	37kDa
Shipping	Ice bag