

## Anti-Pp1 Beta Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62140

Item Number: CM0012829

### Introduction

Rabbit polyclonal antibody against human PP1 beta (PPP1CB) for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse, rat. Essential for cell division, glycogen metabolism, and MAPK signaling.

[Learn More](#)

| Parameter        | Value                                                                                                     |
|------------------|-----------------------------------------------------------------------------------------------------------|
| Product Name     | PP1 beta Rabbit pAb                                                                                       |
| Remarks/Alias    | MP; PP1B; PP1c; NSLH2; PP-1B; PPP1CD; PP1beta; PPP1beta; HEL-S-80p; PP1 beta                              |
| Species          | Human                                                                                                     |
| GeneID (Human)   | 5500                                                                                                      |
| GeneID           | 5500                                                                                                      |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 227-327 of human PP1 beta (NP_002700.1). |
| Source           | Rabbit                                                                                                    |
| Category         | Polyclonal Antibodies                                                                                     |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                  |
| Cross Reactivity | Human, Mouse, Rat                                                                                         |
| SWISS            | P62140                                                                                                    |
| Protein Weight   | 37kDa                                                                                                     |
| Shipping         | Ice bag                                                                                                   |