

Anti-Pp2A-A α /Pr65A/Ppp2R1A Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P30153

Item Number: CM0029476



Introduction

Rabbit polyclonal antibody targeting human PP2A-A α /PR65 α /PPP2R1A, a scaffold subunit of protein phosphatase 2A. Suitable for WB, IF/ICC, and ELISA in Human, Mouse, and Rat samples.

[Learn More](#)

Parameter	Value
Product Name	PP2A-A α /PR65 α /PPP2R1A Rabbit pAb
Remarks/Alias	MRD36; PP2AA; PR65A; PP2AAALPHA; PP2A-Aalpha; PP2A-A α /PR65 α /PPP2R1A
Species	Human
Gene ID (Human)	5518
Gene ID	5518
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human PP2A-Aalpha/PR65 α /PP2A-A α /PR65 α /PPP2R1A (NP_055040.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P30153
Protein Weight	65kDa
Shipping	Ice bag