


Anti-Pp2A-A β /Pr65B/Ppp2R1B Polyclonal Antibody For Wb, Elisa - P30154

Item Number: CM0021567



Introduction

Rabbit polyclonal antibody targeting human PP2A-A β /PR65 β /PPP2R1B, the scaffold subunit of protein phosphatase 2A. Validated for Western blot and ELISA, with cross-reactivity to rat. Suitable for studying PP2A signaling, cell cycle regulation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PP2A-A β /PR65 β /PPP2R1B Rabbit pAb
Remarks/Alias	PPP2R1B; PP2A-Abeta; PR65B; protein phosphatase 2 scaffold subunit Abeta; PP2A-A β /PR65 β /PPP2R1B
Species	Human
Gene ID	5519
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 502-601 of human PP2A-A β /PR65 β /PPP2R1B (NP_002707.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P30154
Protein Weight	52kDa/61kDa/66kDa/73kDa
Shipping	Ice bag