

Anti-Pp2A-B55A/Pr55A/Ppp2R2A Monoclonal Antibody For Wb, Elisa - P63151

Item Number: CM0007200



Introduction

Human PP2A-B55α/PR55α/PPP2R2A rabbit monoclonal antibody validated for WB and ELISA. Recognizes 52kDa protein, cross-reacts with mouse. Ideal for cell cycle and mitosis studies.

[Learn More](#)

Parameter	Value
Product Name	PP2A-B55α/PR55α/PPP2R2A Rabbit mAb
Remarks/Alias	B55A; PR52A; PR55A; B55ALPHA; PR55alpha; PP2A-B55α/PR55α/PPP2R2A
Species	Human
GeneID (Human)	5520
GeneID	5520
Immunogen	Synthetic Peptide - A synthetic peptide corresponding to a sequence within amino acids 94-194 of human PP2A-B55α/PR55α/PPP2R2A (NP_002708.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P63151
Protein Weight	52kDa
Shipping	Ice bag