


Anti-Pp2A-B55A/Pr55A/Ppp2R2A Rabbit Pab - P63151

Item Number: CM0014540



Introduction

High-quality rabbit polyclonal antibody against human PP2A-B55α/PR55α/PPP2R2A (PPP2R2A). Validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PP2A-B55α/PR55α/PPP2R2A Rabbit pAb
Remarks/Alias	B55A; PR52A; PR55A; B55ALPHA; PR55alpha; PP2A-B55α/PR55α/PPP2R2A
Species	Human
GeneID (Human)	5520
GeneID	5520
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-200 of human PP2A-B55α/PR55α/PPP2R2A (NP_002708.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63151
Protein Weight	52kDa
Shipping	Ice bag