


Anti-Pp2A-B55A/Pr55A/Ppp2R2A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P63151

Item Number: CM0020376



Introduction

Rabbit polyclonal antibody targeting PP2A-B55α/PR55α/PPP2R2A, validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Essential for cell cycle regulation and phosphatase research.

[Learn More](#)

Parameter	Value
Product Name	PP2A-B55α/PR55α/PPP2R2A Rabbit pAb
Remarks/Alias	B55A; PR52A; PR55A; B55ALPHA; PR55alpha; PP2A-B55α/PR55α/PPP2R2A
Species	Human
Gene ID (Human)	5520
Gene ID	5520
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-447 of human PP2A-B55α/PR55α/PPP2R2A (NP_002708.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P63151
Protein Weight	52kDa
Shipping	Ice bag