


Anti-Pp2A-B56Δ/Pr61Δ/Ppp2R5D Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q14738

Item Number: CM0003335



Introduction

Rabbit monoclonal antibody targeting human PP2A-B56δ (PPP2R5D). Suitable for WB, IHC-P, ELISA; cross-reacts with Human, Mouse, Rat. Recombinant immunogen (aa 480-602), 70kDa. For research use.

[Learn More](#)

Parameter	Value
Product Name	PP2A-B56δ/PR61δ/PPP2R5D Rabbit mAb
Remarks/Alias	B56D; MRD35; B56delta; PP2A-B56δ/PR61δ/PPP2R5D
Species	Human
GeneID (Human)	5528
GeneID	5528
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 480-602 of human PPP2R5D (NP_006236.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14738
Protein Weight	70kDa
Shipping	Ice bag