


Anti-Pp2A-B56Δ/Pr61Δ/Ppp2R5E Polyclonal Antibody For Wb, Elisa - Q16537

Item Number: CM0020332



Introduction

High-quality rabbit polyclonal antibody targeting human PP2A-B566/PR616 (PPP2R5E) validated for WB and ELISA with cross-reactivity to rat. Ideal for phospho-signaling studies and IVD assay development.

[Learn More](#)

Parameter	Value
Product Name	PP2A-B566/PR616/PPP2R5E Rabbit pAb
Remarks/Alias	B56E; B56epsilon; PP2A-B566/PR616/PPP2R5E
Species	Human
GeneID (Human)	5529
GeneID	5529
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q16537
Protein Weight	55kDa
Shipping	Ice bag