


Anti-Pp2A-B56γ/Pr61γ/Ppp2R5C Rabbit Polyclonal Antibody For Wb, Elisa - Q13362

Item Number: CM0025011



Introduction

Rabbit polyclonal antibody targeting PP2A-B56γ/PR61γ/PPP2R5C, validated for WB and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for studying DNA damage response, ERK signaling, and cell cycle regulation.

[Learn More](#)

Parameter	Value
Product Name	PP2A-B56γ/PR61γ/PPP2R5C Rabbit pAb
Remarks/Alias	B56G; PR61G; B56gamma; PP2A-B56γ/PR61γ/PPP2R5C
Species	Human
GeneID (Human)	5527
GeneID	5527
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-320 of human PP2A-B56γ/PR61γ/PPP2R5C (NP_848702.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13362
Protein Weight	61kDa
Shipping	Ice bag