

Anti-Pp2A-B56B/Pr61B/Ppp2R5B Polyclonal Antibody For Wb And Elisa - Q15173

Item Number: CM0021558

Introduction

Rabbit polyclonal antibody targeting human PP2A-B56β/PR61β/PPP2R5B, validated for WB and ELISA with rat cross-reactivity. Ideal for PP2A regulation, AKT1 signaling, and cytoplasmic phosphatase studies.

[Learn More](#)

Parameter	Value
Product Name	PP2A-B56β/PR61β/PPP2R5B Rabbit pAb
Remarks/Alias	B56B; PR61B; B56beta; PP2A-B56β/PR61β/PPP2R5B
Species	Human
GeneID (Human)	5526
GeneID	5526
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-497 of human PP2A-B56β/PR61β/PPP2R5B (NP_006235.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	Q15173
Protein Weight	57kDa
Shipping	Ice bag