


Anti-Pp2A Catalytic A Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc - P67775

Item Number: CM0028885



Introduction

PP2A Catalytic α rabbit pAb (P67775) validated for WB, IHC-P, IF/ICC, IP, ELISA with human, mouse, rat cross-reactivity. Ideal for studying PP2A signaling, cell cycle regulation, and transcription control.

[Learn More](#)

Parameter	Value
Product Name	PP2A Catalytic α Rabbit pAb
Remarks/Alias	RP-C; PP2Ac; PP2CA; NEDLBA; PP2Calpha; PP2A Catalytic α
Species	Human
GeneID (Human)	5515
GeneID	5515
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-309 of human PP2A Catalytic α (NP_002706.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P67775
Protein Weight	36kDa
Shipping	Ice bag