

Anti-Pp2A Catalytic B (Ppp2Cb) Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P62714

Item Number: CM0025020

Introduction

Rabbit polyclonal antibody targeting human PP2A catalytic subunit beta isoform (PPP2CB). Validated for WB, IHC-P, IF/ICC, ELISA; cross-reacts with mouse and rat. Ideal for studying PP2A signaling and STRIPAK complexes.

[Learn More](#)

Parameter	Value
Product Name	PP2A Catalytic β Rabbit pAb
Remarks/Alias	PP2CB; PP2Abeta; PP2A Catalytic β
Species	Human
GenelD (Human)	5516
GenelD	5516
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-309 of human PP2A Catalytic β (NP_001009552.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P62714
Protein Weight	36kDa
Shipping	Ice bag