


Anti-Ppa1 Polyclonal Antibody For Wb, Elisa - Q15181

Item Number: CM0018362



Introduction

High-quality rabbit polyclonal antibody against human PPA1 (Inorganic pyrophosphatase). Validated for Western blot and ELISA; cross-reacts with mouse and rat. Ideal for studying phosphate metabolism and JNK signaling pathways in research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	PPA1 Rabbit pAb
Remarks/Alias	PP; PP1; IOPPP; HEL-S-66p; SID6-8061; PPA1
Species	Human
GeneID (Human)	5464
GeneID	5464
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-289 of human PPA1 (NP_066952.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15181
Protein Weight	33kDa
Shipping	Ice bag